


Face mask and respirator materials list

➤ Surgical masks

Product Picture



Material

A	Round Elastic Band
B	Plastic Nose Wire
1	Blue Non-woven Fabric
2	Filtering Material (Melt-Blown)
3	White Non-woven Fabric

➤ Dust Masks

Product Picture





Material

A	Flat Elastic Band
B	Plastic Metal Nose Wire
1	PP Non-woven Fabric
2	Filtering Material (Melt-Blown)
3	Active Carbon Non-woven
4	PP Non-woven Fabric

➤ Respirator

Product Picture



Material

A	Headband Elastic
B	Aluminum Nose Clip
C	Nose Foam
D	Mask Ear-loop Holder
1	Non-woven Fabric
2	Filtering Material (Melt-Blown)
3	Forming Cotton

Other Material

Exhalation Valve

